

⚠ DANGER

Always disconnect the power cord from the wall outlet before performing service on the table to prevent electrical shock. Failure to do so could result in severe personal injury or death.

⚠ WARNING

Refer to the Operator's Manual supplied with the unit for complete instructions on operation. Failure to do so could result in personal injury.



Qty.	Description
1	Decorative Side Shroud
1	Silicone Bumper

The Brewer Company, LLC
N88 W13901 Main Street, Suite 100
Menomonee Falls, WI 53051
P 1.888.BREWER.1
F 1.262.251.2332
www.brewercompany.com

Tools Required

1.1/8" Allen Wrench

Note: The decorative side shroud on your model may look different from the one pictured below; however, the procedure is still the same.

Decorative Side Shroud Removal:

1. Raise seat to the full up position.
2. Tilt seat to the full back position.

Note: These above two steps are not necessary for changing a decorative side shroud, but will make the procedure much easier.


3. Unplug the table from wall outlet.
4. If present, remove clear silicone bumper at top of decorative side shroud.



5. Grasp decorative side shroud with both hands and lift straight up. Pull decorative side shroud toward you as you "wiggle" shroud slightly to allow screw heads to pass through the (4) keyhole slots that secure decorative side shroud in place.



6. Note the (4) keyhole slots into which the new shroud screw heads will pass.



Decorative Side Shroud Installation:

1. Align the four screw heads with the four keyhole slots. Insert the screw heads into the keyhole slots and gently push down. Note: To avoid damage do not force the decorative side shroud into place. If decorative side shroud does not attach easily go to step 2, otherwise proceed to step 3



2. The decorative side shroud comes from the factory preadjusted and should fit, or be very close. If fit is too tight loosen the allen screw on the corner that is tight by 1/8 turn. Try fit again. Continue this procedure until it fits correctly.



3. Replace the clear silicone bumper at the top of the decorative side shroud to ensure that it holds the shroud in place. (Make sure bumper is snug to shroud to prevent shroud from popping off under normal use)
4. Plug table back into wall outlet.
5. Using foot control or hand pendant tilt seat forward until it stops.
6. Lower and raise table two times to ensure that the shrouds are functioning properly.

