


⚠ DANGER

Always disconnect the power cord from the wall outlet before performing service on the table to prevent electrical shock. Failure to do so could result in severe personal injury or death.

⚠ WARNING

Refer to the Operator's Manual supplied with the unit for complete instructions on operation. Failure to do so could result in personal injury.



Description	Qty.
Foot Control	1

The Brewer Company
13901 Main Street
Menomonee Falls, WI 53051
P 1.888.BREWER.1
F 1.262.251.1786
www.brewercompany.com

INSTALLATION

1. Unplug the old Foot Control from the receptacle at the back of the Power Table's base.
2. Plug the replacement Foot Control into the receptacle aligning the keyway. Make sure the plug is fully inserted into the receptacle.

